### **Product Features**

4 independent channels with 8 isolated outputs

Laser current source modules from 200mA to 8A

LD controller modules from 200mA to 2A with integrated 12W TEC

32W TEC only modules with voltage measurement

TEC modules accept thermistor, IC and RTD temperature sensors

Standard GPIB IEEE interface

The LDC-3900 Modular Laser Diode Controller features four channels with eight isolated outputs for controlling multiple laser diodes. Modules include current sources with maximum outputs from 200mA to 8A, TE (thermoelectric) controller modules with up to 32W of power and voltage measurement, and controller modules with laser control from 200mA to 2A and integrated 12W TE control.

An independent power supply powers each channel, providing protection and stability for your laser diode. Every module incorporates low noise, high stability performance, and ILX Lightwave's unmatched laser diode protection topologies including independent current limits and laser diode shorting relays.

Remote communication through an IEEE/GPIB interface simplifies testing and control of multiple devices. LabVIEW® drivers are also available for any combination of mainframe and module.



Modular Laser Diode Controller



Flexible, Comprehensive Control of Laser Diodes



# LDC 3900

Modular Laser Diode Controller

## Wide Range of Modules

Five current source modules and five combination modules along with two TE controller only modules make the LDC-3900 configurable for many laser diode testing and control applications. Each module is electrically floating or fully isolated from all other modules. This allows you to configure your laser diode test system without the worry of potential laser diode damaging ground loops.

#### **Current Source Modules**

The LDC-3900 current source design offers superior laser protection and low noise, high stability performance. These modules also feature a photo-diode measurement circuit for devices with backfacet monitor photodiode and analog modulation up to 500kHz for dithering the laser current for wavelength tuning. Five different current source modules up to 8A can be driven in any one of the following modes:

- 1) Constant current, low bandwidth
- 2) Constant current, high bandwidth
- 3) Constant optical power

## **Highly Stable Temperature Control**

The LDC-3900 TEC modules control temperature of your devices with 32W of power. These TEC modules offer maximum flexibility with a choice of operating modes and temperature sensors covering thermistors, IC, and RTDs. A low noise, biopolar output with TE voltage measurement and an ultra-stable topology achieves stabilities better than 0.005°C. A smart integrator control loop, programmable from the front panel or through GPIB, delivers fast settling times.

### **Controller Modules**

Controller modules combine a current source with a temperature controller in one module. Laser current from 200mA to 2A is available with an integrated 12W temperature controller for current and temperature control of laser diodes. All of the features found in the current only and temperature control modules are incorporated into these flexible modules and include multiple modes of operation, external modulation, a choice of temperature sensors, and all protection features.

### **Intuitive Front Panel**

Divided into two sections, TEC and LASER, the front panel offers quick, easy operation and information display without confusing multi-layer menus. Each channel is directly addressable from the front panel "adjust" section and indicated through discrete LEDs next to the respective display. Laser and TEC parameters and modes are easily selected or adjusted through discrete push buttons.

#### Powerful GPIB Interface

For automated control, the IEEE/GPIB interface allows programming and readout from most computers. All instrument and module functions are accessible through the interface allowing you to simultaneously control multiple laser diodes from the same address. For virtual instrument programming, LabVIEW® drivers are available upon request or through the ILX website.

\* Semiconductor lasers are sensitive devices. Always take appropriate antistatic precautions and use extreme care when handling laser diodes. For more information, request ILX Application Note #3, "Protecting Your Laser Diode."



Up to four modules can be easily adjusted and controlled from the LDC-3900 front panel.

## **Specifications**

## Current Source Modules 1

ourrent bource mou	uico				
Current Source 1	39020	39050	39100	39400	39800
DRIVE CURRENT OUTPUT Output Current Range: Setpoint Resolution: Setpoint Accuracy: Compliance Voltage:	0–200 mA 10 μA ±0.1% of FS 7V	0–500 mA 10 μA ±0.1% of FS 6.5V	0–1000 mA 100 μA ±0.1% of FS 6V	0–4000 mA 100 μA ±0.1% of FS 5V	0–8000 mA 125 μA ±0.1 % of FS 5V at connector
Temperature Coefficient: Short-Term Stability (one hour): <sup>2</sup> Long-Term Stability (24 hours): <sup>3</sup> Noise and Ripple (μA rms) <sup>4</sup>	<60 ppm/°C <20 ppm <50 ppm	<60 ppm/°C <20 ppm <40 ppm	<100 ppm/°C <20 ppm <40 ppm	<100 ppm/°C <20 ppm <40 ppm	(4.5V at end of cable) <100 ppm/°C <20 ppm <40 ppm
High Bandwidth Mode: Low Bandwidth Mode: With LNF-320: <sup>5</sup> Transients:	<4 μA <3 μA <1 μA	<5 μA <3 μA <2 μA	<10 μA <5 μA <2.5 μA	<20 μA <20 μA <5 μA	<120 μA <110 μA N/A
Operational: <sup>6</sup> Power-line spike induced: <sup>7</sup>	<1 mA <5 mA/<8 mA	<1 mA <5 mA/<8 mA	<2 mA <5 mA/<8 mA	<5 mA <10 mA/<20 mA	<8 mA <20 mA/<40 mA
Isolation:		All modules isolated	from other modules ar	nd earth ground	
DRIVE CURRENT LIMIT SETTI Range: Resolution: Accuracy:	NGS 0–200 mA 0.5 mA ±2 mA	0–500 mA 2 mA ±5 mA	0–1000 mA 4 mA ±10 mA	0–4000 mA 16 mA ±40 mA	0-8000 mA 40 mA ±80 mA
PHOTODIODE FEEDBACK Type: Reverse Bias: Photodiode Current Range: Output Stability: 8 Setpoint Accuracy:	Transimpedance 0-5V, adjustable 0-5 mA 0.02% ±0.05% of FS	Transimpedance 0-5V, adjustable 0-5 mA 0.02% ±0.05% of FS	Transimpedance 0–5V, adjustable 0–10 mA 0.02% ±0.05% of FS	Transimpedance 0-5V, adjustable 0-20 mA 0.02% ±0.1% of FS	Transimpedance 0-5V, adjustable 0-20 mA 0.02% ±0.1% of FS
EXTERNAL ANALOG MODULAT Input: Transfer Function: Bandwidth (3dB) High Bandwidth: Low Bandwidth: Low Bandwidth:	TON 0–10V, 10kΩ 20 mA/V DC to 500kHz DC to 5KHz DC to 30Hz	0–10V, 10kΩ 50 mA/V DC to 200kHz DC to 5kHz DC to 30Hz	0–10V, 10kΩ 100 mA/V DC to 200kHz DC to 5kHz DC to 30Hz	0–10V, 10kΩ 400 mA/V DC to 50kHz DC to 2kHz DC to 30Hz	0–10V, 10kΩ 800 mA/V DC to 50kHz DC to 2kHz N/A
OUTPUT CONNECTORS Current Source Output: Photodiode Input: Analog Modulation Input:	9-pin, D-sub Coax BNC Coax BNC inst. amp. input	9-pin, D-sub Coax BNC Coax BNC inst. amp. input	9-pin, D-sub Coax BNC Coax BNC inst. amp. input	9-pin, D-sub Coax BNC Coax BNC inst. amp. input	16-pin, Bulkhead Coax BNC Coax BNC inst. amp. input
MEASUREMENT (DISPLAY) 10 Output Current Range: Resolution: Accuracy: 11	0–200 mA 0.01 mA ±0.05% of FS	0-500 mA 0.01 mA ±0.1% of FS	0-1000 mA 0.1 mA ±0.1% of FS	0-4000 mA 0.1 mA ±0.1% of FS	0-8000 mA 0.1 mA ±0.1% FS
Photodiode Current Range: Resolution: Accuracy: 11 Photodiode Proposition	0–5000 μA 1 μA ±0.05% of FS	0–5000 μA 1 μA ±0.05% of FS	0–10,000 μA 1 μA ±0.05% of FS	0–20,000 μA 1 μA ±0.05% of FS	0–20,000 μA 1 μA ±0.05% of FS
Photodiode Responsivity Range: 12 Resolution: Optical Power	0.01 μA/mW ·	0.01 μA/mW ·	0.00–600.00 μA/mW 0.01 μA/mW	0.01 μA/mW ·	0.00–1000.00 μA/mW 0.01 μA/mW
Range: Resolution: Forward Voltage	0.00–200.00 mW 0.01 mW	0.00–500.00 mW 0.1 mW	0.00–1000.0 mW 0.1 mW	0.00–5000.0 mW 0.1 mW	0.00–8000.0 mW 0.1 mW
Range: Resolution: Accuracy: <sup>13</sup>	0.000-7.000V 1 mV ±5 mV	0.000-7.000V 1 mV ±5 mV	0.000- 7.000V 1 mV ±5 mV	0.000–5.000V 1 mV ±5 mV	0.000-5.000V 1 mV ±5 mV

#### **CURRENT SOURCES NOTES**

- URRENT SOURCES NOTES
  All values relate to a one-hour warm-up period.
  Over any non-hour period, half-scale output at 25°C ambient.
  Over any 24-hour period, half-scale output at 25°C ambient.
  Measured optically from resulting intensity fluctuations of a laser diode with a 150kHz bandwidth photodetector. Measurements made with 1MHz detector are typically 10% higher.
  With model LNF-320 low noise CW filter enabled.
  Maximum output current transient resulting from normal operational situations (i.e., power ine plug removal). For more information, request ILX "Transient Test Standards" #LDC-00196.
- Maximum output current transient resulting from a 1000V power line transient spike.
   Tested to ILX Lightwave Technical Standard #LDC-00196.
   Maximum monitor photodiode current drift over any 30 minute period. Assumes

- waxman more processes and the second content and the second con

12 Responsivity value is user-defined and is used to calculate the optical power.

13 Voltage measurement accuracy while driving calibration load. Connected at the rear panel connector. Accuracy may vary depending on load and cable length used.

Modular Laser Diode Controller

Modular Laser Diode Controller

### **Specifications**

#### TEC Modules <sup>1</sup>

#### TCM-39034 15

TEMPERATURE CONTROL

Temperature Control Range: 2 -99.9°C to 199.9°C Thermistor Setpoint Resolution Resolution and Accuracy: 3 Accuracy 3 -20°C to 20°C 0.1°C ±0.2°C 20°C to 50°C 0.2°C ±0.2°C

AD590 and LM335 Setpoint Resolution and Accuracy:

Resolution Accuracy 3 0.01°C -20°C to 50°C ±0.2°C Short Term Stability (one hour.): 5 Long Term Stability (24 hours): 6 <±0.004°C <±0.01°C

TEC OUTPUT 7

Output Type: Bipolar, constant current source Isolation: Isolated from other modules and earth ground

Compliance Voltage: >8V Short Circuit Output Current: 4A Maximum Output Power: 32W Current Limit Range: 0-4A Current Limit Set Accuracy: ±50mA Ripple/ Noise: 8 <1mA, rms

Smart Integrator, Hybrid PI Control Algorithm:

#### TEMPERATURE SENSOR

Thermistor: IC Temperature Sensor: RTD Sensor: 9

Thermistor Sensing Current: Sensor Bias:

Usable Thermistor Range: Typical Sensor Output 1 AD590 Current Output: LM335 Voltage Output: RTD (PT100) Resistance: User Calibration:

NTC (2-wire) AD590/LM335 Pt100/other 100Ω RTD 10/100uA

AD590 = 8V, LM335 = 1mA, RTD = 0.8 mA 9  $25-450,000\Omega$ , typical

 $\begin{array}{l} I~(25^{\circ}C) = 298.2 \mu A, \ It = 1 \mu A/K \\ V~(25^{\circ}C) = 2.73 V, \ Vt = 10 m V/K \end{array}$ R (25°C) =  $109.73\Omega$ Thermistor = Steinhart-Hart IC Sensors, RTD = Two-point

#### TEC OUTPUT CONNECTORS

Temperature Controller Output: 15-pin, D-sub

#### TEC MEASUREMENT (DISPLAY) 11

10µA Setting: 13 100µA Setting: 14 0.01°C -99.99°C to 199.99°C ±0.05°C Thermistor Resistance: 10µA Setting: ±0.05% 0.0-480.00k $\Omega$  $0.01k\Omega$ 

-99.99°C to 199.99°C

Accuracy

+0.1°C

Resolution

0.01°C

100µA Setting:  $0.0 - 48.000 k\Omega$  $0.001 k\Omega$ ±0.005% TE Current: -4.000 to 4.000A 0.001A ±0.04A

TE VOLTAGE 15

Temperature:

Voltage Range: -9.999 to 9.999V Voltage Resolution: 1mV ±30mV 16 Voltage Accuracy:

#### TEC MODULE NOTES:

All values relate to a one-hour warm-up period.

Software limits of range. Actual range possible depends on the physical load, thermistor type, and TE module used.

thermistor type, and 1 E module used.

Accuracy figures are quoted for a typical 10kΩ thermistor and 100µA current setting. 
Accuracy figures are relative to the calibration standard. Both resolution and accuracy 
are dependent upon the user-defined configuration of the instrument. 
Accuracy depends upon the sensor model selected, the calibration standard,

Accuracy depends upon the sension indeed selected, the calibration standard, and the user-defined configuration of the instrument. Over any one-hour period, half-scale output, controlling an LDM-4412 mount at 25°C, with  $10 \mathrm{k}\Omega$  thermistor on  $100 \mathrm{\mu}\Delta$  setting. Over any 24-hour period, half-scale output, controlling an LDM-4412 mount at 25°C, with  $10 \mathrm{k}\Omega$  thermistor on  $100 \mathrm{\mu}\Delta$  setting.

6

9

with 10x2 definision of 100pA seams.

Into a 152 load.

Measured at 1A output over a bandwidth of 10Hz to 10MHz.

To use RTD sensors with model 39034, order TSC-599 Temperature Sensor Converter accessory.

Nominal temperature coefficients, It and Vt, apply over the rated IC temperature

10

sensor range.
Displayed on LDC-3900 mainframe front panel "TEC" section. 11

Displayed on LDC-3900 maintraine front panel TEC section. Software limits of display range. Using a 100k $\Omega$  thermistor, controlling an LDM-4412 mount over ~30°C to 25°C. Using a 10k $\Omega$  thermistor, controlling an LDM-4412 mount over 0°C to 90°C. Model 39034 has TEC Voltage measurement through GPIB only. Voltage measurement accuracy while driving calibration load. Accuracy is dependent upon load used.

#### MAINFRAME/GENERAL

Chassis Ground: 4mm banana jack
Power Requirements, 50–60Hz: 100VAC, 120VAC, 220VAC, 240VAC, (+6%/-10%)

(user selectable) 145mm x 426mm x 346mm Size (HxWxD): (5 5/8" x 16 3/4" x 13 5/8") 12.5kg (27.5lbs) 16.5 kg (37.7lbs), typical 0°C to 40°C Weight Mainframe: Weight Mainframe (loaded): Operating Temperature: Storage Temperature: -40 to 70°C

<90%, non-condensing Keyswitch, Interlock, Output Delay Laser Safety Features: (meets CDRH US21 CFR 1040.10)

Laser Display: 5-digit, green LED TEC Display: 5-digit, green LED

#### ORDERING INFORMATION NOTES

Includes ILX model LTS-520 calibrated 10kΩ thermistor

In keeping with our commitment to continuing improvement, ILX Lightwave reserves the right to change specifications at any time without notice and with out liability for such changes.

#### ORDERING INFORMATION

LDC-3900 Modular Laser Diode Controller Mainframe 200mA Current Source Module CSM-39020 500mA Current Source Module 1A Current Source Module CSM-39050 CSM-39100 CSM-39400 4A Current Source Module

CMS-39800 8A Current Source Module (Module take two slots in LDC-3900) TCM-39034\*

32W TEC Module with Voltage Measurement
Current/TEC Combination Module (200mA Drive Current/8W TEC)
Current/TEC Combination Module (500mA Drive Current/12W TEC) I CM-39420 LCM-39425 Current/TEC Combination Module (500mA Drive Current/12W TEC) with Modulation LCM-39427

LCM-39437 Current/TEC Combination Module (1A Drive Current/12 W TEC) with

Modulation I CM-39440 Current/TEC Combination Module (2A Drive Current/8W TEC)
Current Source/Laser Diode Mount Interconnect Cable CC-305S CC-306S CC-501S Current Source/Unterminated Interconnect Cable TE Controller/Unterminated Interconnect Cable CC-505S TE Controller/Laser Diode Mount Interconnect Cable

Calibrated 10kΩ Thermistor Uncalibrated 10kΩ Thermistor TS-520 TS-530

Uncalibrated AD590LH IC Temperature Sensor Uncalibrated AD590LH IC Temperature Sensor Uncalibrated LM355AH IC Temperature Sensor RTD Temperature Sensor Control Option (for 39032 Module) Rack Mounting Kit TS-540 TSC-595

RM-103 Unipolar Heater Control Adapter Low Noise Filter LICA-350

LNF-320 LabVIEW®Instrument Driver



31950 Frontage Road, Bozeman, MT 59715 • FAX: 406-586-9405

www.newport.com/ilxlightwave





## Specifications

## **Combination Modules**

	20/20	20/25	20/2=	20/2=	20//0	
	39420	39425	39427	39437	39440	
ISOLATION	Each module is isolated from other modules and earth ground. TEC and current source independently isolated					
OUTPUT CONNECTORS Current Source Output: Photodiode Input:	9-pin, D-sub 15-pin, D-sub	9-pin, D-sub 15-pin, D-sub	9-pin, D-sub 15-pin, D-sub	9-pin, D-sub 15-pin, D-sub	9-pin, D-sub 15-pin, D-sub	
DRIVE CURRENT OUTPUT Output Current Range: 9 Setpoint Resolution: Setpoint Accuracy: Compliance Voltage:	0–200 mA 10 μA ±0.05% of FS 5V	0–500 mA 10 μA ±0.05% of FS 6V	0–500 mA 10 μA ±0.05% of FS 6V	0–1000 mA 100 μA ±0.01% of FS 6V	0–2000 mA 100 μA ±0.01% of FS 5V	
Temperature Coefficient: Short-Term Stability (one hour): <sup>2</sup> Long-Term Stability (24 hours): <sup>3</sup> Noise and Ripple (μA/rms) <sup>4</sup>	100 ppm/°C 25 ppm 50 ppm	100 ppm/°C 25 ppm 50 ppm	100 ppm/°C 25 ppm 50 ppm	100 ppm/°C 25 ppm 50 ppm	100 ppm/°C 25 ppm 50 ppm	
Unfiltered: With model 320 Filter: 5 Transients:	<2.5 μA <1 μA	<4 μA <1.5 μA	<4 μA <1.5 μA	<4 μA <1.5 μA	<10 μA <2 μA	
Operational: 6 1kV EFT: Surge: 7  DRIVE CURRENT LIMIT S	<1 mA <4 mA <7 mA	<1 mA <4 mA <7 mA	<1 mA <4 mA <7 mA	<1 mA <4 mA <7 mA	<1 mA <10 mA <8 mA	
Range: Accuracy:	0–200 mA ±2 mA	0–500 mA ±5 mA	0–500 mA ±5 mA	0–1000 mA ±10 mA	0–2000 mA ±20 mA	
PHOTODIODE FEEDBACK Type: Range: Output Stability: 8 Setpoint Accuracy:	20–2000 μA ±2 μA ±2 μA	Current input differer 20–2000 µA ±2 µA ±2 µA	ntial, zero bias, all modu 20–4000 μΑ N/A N/A	iles 20–4000 μΑ N/A N/A	20–5000 μA ±2 μA ±5 μA	
EXTERNAL ANALOG MODI Input: Transfer Function: Bandwidth (3dB):	ULATION N/A N/A N/A	N/A N/A N/A	0–10V, 10kΩ 50 mA/V DC to 250 kHz	0–10V, 10kΩ 100 mA/V DC to 200 kHz	N/A N/A N/A	
DRIVE CURRENT MEASUF Output Current Range: Output Current Resolution: Output Current Accuracy: 10 Photodiode Current Range: PD Current Resolution: PD Responsivity Range: PD Responsivity Resolution: Optical Power Range: Optical Power Resolution: Voltage Range: Voltage Accuracy:	0-200.00 mA 0.01 mA ±0.1 mA ±0.1 mA 0-2000 μA 1 μA 0.00-1000.00μA/mW 0.01 μA/mW 0.00-200.00 mW 10 μW N/A N/A	0–500.00 mA 0.01 mA ±0.5 mA 0–2000 μA 1 μA	0–500.00 mA 0.01 mA ±0.5 mA 0–4000 μA 1 μA / 0.00–1000.00μA/mW 0.01 μA/mW 0.00–1000.00 mW 10 μW 0 to 7V ±5mV	0-1000.00 mA 0.01 mA ±1 mA 0-4000 μA 1 μA 0.00-1000.00μA/mW 0.01 μA/mW 0.00-1000.00 mW 10 μW 0 to 7V ±5mV	0- 2000.0 mA 0.1 mA ±1 mA 0-5000 μA 1 μA 0.00-1000.00μA/mW 0.01 μA/mW 0.00-2000.0 mW 100 μW N/A	
TEMPERATURE CONTROL Temperature Control Range: 11 Thermistor Setpoint	–99.9°C to 99.9°C	-99.9°C to 99.9°C	−99.9°C to 99.9°C	−99.9°C to 99.9°C	−99.9°C to 99.9°C	
Resolution and Accuracy <sup>12</sup> -20°C to +20°C: +20°C to +50°C: Short Term Stability (1 hr.): <sup>13</sup> Long Term Stability (24 hrs.): <sup>14</sup>	Res. Acc. 0.1°C ±0.2°C 0.2°C ±0.2° C <±0.05°C <±0.1° C	Res. Acc. 0.1°C ±0.2°C 0.2°C ±0.2°C <±0.05°C <±0.1°C	Res. Acc. 0.1°C ±0.2°C 0.2°C ±0.2°C <±0.05°C <±0.1°C	Res. Acc. 0.1°C ±0.2°C 0.2°C ±0.2°C <±0.05°C <±0.1°C	Res. Acc. 0.1°C ±0.2°C 0.2°C ±0.2°C <±0.05°C <±0.1°C	
Output Type: Compliance Voltage: Short Circuit Output Current: Maximum Output Power: Current Noise and Ripple: Current Limit Range: Current Limit Set Accuracy: Control Algorithm:	>4V DC 2 A 8 W <1 mA rms 0–2 A 0.05 A	Bipolar, constant cur >6V DC 2 A 12 W <1 mA rms 0-2 A 0.05 A Smart Integrator, Hyl	rent source, all modules >6V DC 2 A 12 W <1 mA rms 0-2 A 0.05 A brid PI, all modules	5 >6V DC 2 A 12 W <1 mA rms 0-2 A 0.05 A	>4V DC 2 A 8 W <1 mA rms 0–2 A 0.05 A	

# LDC 3900

## Modular Laser Diode Controller

# LDC 3900

## Modular Laser Diode Controller

## **Specifications**

	39420	39425	39427	39437	39440	
TEMPERATURE SENSOR						
Types:	Thermistor (2-wire NTC)					
Thermistor Sensing Current:	10/100 µA (user-selectable)	10/100 µA (user-selectable)	10/100 μA (user-selectable)	10/100 μA (user-selectable)	10/100 μA (user-selectable)	
Usable Thermistor Range:	25–450,000 $\Omega$ typical					
User Calibration:	Steinhart-Hart	Steinhart-Hart	Steinhart-Hart	Steinhart-Hart	Steinhart-Hart	
TEC MEASUREMENT (DISPLAY)						
Range						
Temperature:	−99.9°C to 99.9°C	−99.9 °C to 99.9°C	−99.9 °C to 99.9°C	−99.9 °C to 99.9°C	−99.9 °C to 99.9°C	
Thermistor Resistance						
10 μA Setting:	0.00–450.00 kΩ	$0.00-450.00 \text{ k}\Omega$	0.00–450.00 kΩ	0.00–450.00 kΩ	0.00–450.00 kΩ	
100 μA Setting:	$0.000$ – $45.000$ k $\Omega$	$0.000-45.000 \text{ k}\Omega$	$0.000-45.000 \text{ k}\Omega$	0.000-45.000 kΩ	0.000–45.000 kΩ	
TE Current:	-2.000 to 2.000 A					
Accuracy						
Temperature:	±0.5°C	±0.5°C	±0.5°C	±0.5°C	±0.5°C	
Thermistor Resistance						
10 μA Setting:	$\pm 0.05  \mathrm{k}\Omega$	$\pm 0.05~\text{k}\Omega$	$\pm 0.05 \text{ k}\Omega$	$\pm 0.05~\text{k}\Omega$	$\pm 0.05~\text{k}\Omega$	
100 μA Setting:	$\pm 0.005~\text{k}\Omega$	$\pm 0.005~\text{k}\Omega$	$\pm 0.005~\text{k}\Omega$	$\pm 0.005~\text{k}\Omega$	$\pm 0.00~5 \mathrm{k}\Omega$	
TE Current:						

#### **COMBINATION MODULES NOTES:**

- 1 All values measured after a one-hour warm-up period.
- 2 Over any one hour period, half-scale output.
- 3 Over any 24 hour period, half-scale output.
- 4 Measured from resulting intensity fluctuations of a laser diode, measured optically with a 150 kHz bandwidth photo detector. Measurements made with 1 MHz detector are typically 10% higher.
- 5 ILX Lightwave model LNF-320 low-noise filter option may be used if lower noise performance is required.
- 6 Maximum output current transient resulting from normal operational situations (e.g., power on-off), as well as accidental situations (e.g., power line plug removal). For more information request ILX "Transient Test Standards" #LDC-00196.
- 7 Maximum output current transient resulting from a 1000V power line transient spike. Tested to ILX Lightwave Technical Standard #LDC-00196.
- 8 Maximum monitor photodiode current drift over any 30-minute period. Constant-power mode stability specification assumes zero drift in detector responsivity.
- 9 Output current rated into a  $1\Omega$  load.
- 10 Measured at 25°C.
- 11 Software limits of range.
- 12 Accuracy figures quoted for a 10kΩ thermistor. Accuracy figures are relative to calibration standard. Both resolution and accuracy are dependent on the user-defined configuration of the instrument.
- 13 Over any one-hour period at 25°C. Short-term temperature stability is a strong function of the thermal environment of the thermistor and TE module. Room air currents in particular can easily cause fluctuations of 0.1°C in an exposed mounting configuration.
- 14 Over any 24 hour period at 25°C.

In keeping with our commitment to continuing improvement, ILX Lightwave reserves the right to change specifications without notice and without liability for such changes.





