


# Connecting the TRA-V6 in a vacuum chamber to a Newport controller



Although the TRA6PPV6, TRA12PPV6 and TRA25PPV6 actuators are vacuum compatible, the SUB-D25M connector supplied with the actuators is **not** vacuum compatible or intended for use in a vacuum environment. Therefore, **customers are responsible** for cutting the supplied cable into 2 sections and making connections to a vacuum feedthrough/bulkhead coupling to connect the TRA-V6 actuator to a Newport controller. This note describes the steps for this procedure.

1. Cut the supplied cable into 2 sections.

The cable consists of 2 bundles, “T1” and “T2”. T1 contains 7 wires and T2 contains 4 wires, as shown in this cross-section:



Each wire's bundle, color, function, wire gauge and pin # on the SUB-D25M connector is described in this table:

| Bundle | Color  | Function       | AWG | SUB-D25M Pin # |
|--------|--------|----------------|-----|----------------|
| T1     | Green  | Ground         | 27  | 22             |
| T1     | Brown  | +5V Supply     | 27  | 21             |
| T1     | Black  | [not used]     | 27  | 20             |
| T1     | Violet | [not used]     | 27  | 19             |
| T1     | Grey   | Negative Limit | 27  | 18             |
| T1     | Pink   | Positive Limit | 27  | 17             |
| T1     | Blue   | Origin         | 27  | 13             |
| T2     | White  | + Phase 2      | 26  | 5              |
| T2     | Orange | – Phase 2      | 26  | 7              |
| T2     | Red    | + Phase 1      | 26  | 1              |
| T2     | Green  | – Phase 1      | 26  | 3              |

**Note: Both T1 and T2 bundles contain a green wire. Be sure to identify and distinguish each wire for appropriate connections.**

2. Connect the flying leads of the cable section attached to the TRA-V6 (inside the vacuum chamber) to the feedthrough.
3. Connect the flying leads of the cable section attached to the SUB-D25M connector to the feedthrough portion outside the vacuum chamber. **IMPORTANT:** Keep the SUB-D25M connector attached to the cable. All EEPROM data is stored in the connector, allowing for quick detection of the stage without manual configuration.
4. Connect the SUB-D25M connector to the appropriate Newport controller.

For more information about the TRA Series Actuators, please see the [TRA User's Manual](#).