



MODEL: LPXY-1H OR M-LPXY-1H

MODEL	D in (mm)	MOUNTING HOLE THREAD
LPXY-1H	0.50 (12.7)	8-32
M-LPXY-1H	0.47 (12.0)	M4 x 0.7

TECH SPECS:

1. X / Y TRANSLATION: ±1.6mm AT 254µm/rev; 0.71µm SENSITIVITY
2. CLEAR APERTURE: Ø.90[22.9mm]
3. OPTICS AXIS HEIGHT: 1.0[25.4mm]
4. MAX OPTICS EDGE THICKNESS: 0.40 [10.2mm]

