



ACCEPTS DIA .375 (9.5)
ACTUATORS

NOTES:

1. TRACK LOCK SCREW AND MICROMETER CLAMP ACCEPT 5/64" (2.0) ALLEN WRENCH
2. ENGLISH MOUNTING HOLES ON 1/2" GRID. METRIC MOUNTING HOLES ON 12.5MM GRID
3. FOR MOUNTING HOLE LOCATION REFER TO DRAWING 9065-X

1.00 (25.4)

2X ALIGNMENT HOLES FOR
5/64 (2.0) ALLEN WRENCH

1/4-20 (M6) THD

8X 8-32 (M4) THD

1.55 (39.4)

.78 (19.7)

.78 (19.7)

1.55 (39.4)

.60 (15.1)

2.45 (62.1)

9065-Z(M) TRANSLATION STAGE

 New Focus™