




77565 BIFURCATED H.G. FUSED SILICA BUNDLE E



77533 BIFURCATED GLASS BUNDL E



77536 TRIFURCATED GLASS BUNDLE



Fiber Material	N.A.	Accept. Cone 2a (degrees)	Common Bundle Dia. in (mm) (A)	Common End Tip Dia. in (mm) (B)	Common End Tip Length in (mm) (C)	Total Length in (mm) (D)	Leg Bundle Dia. in (mm) (E)	Leg End Tip Dia. in (mm) (F)	Leg End Tip Length in (mm) (G)	Model No.
<b>Bifurcated Bundles</b>										
Glass	0.56	68	0.177 (4.5)	0.250 (6.4)	0.50 (12.7)	36 (914)	0.125 (3.2)	0.187 (4.75)	0.50 (12.7)	<b>77533</b>
H.G. Fused Silica	0.22	25	0.125 (3.2)	0.187 (4.7)	0.50 (12.7)	36 (914)	0.088 (2.2)	0.187 (4.75)	0.50 (12.7)	<b>77565</b>
<b>Trifurcated Bundles</b>										
Glass	0.56	68	0.281 (5.5)	0.312 (7.9)	0.50 (12.7)	36 (914)	0.125 (3.2)	0.187 (4.75)	0.50 (12.7)	<b>77536</b>