



GLASS AND STANDARD GRADE FUSED SILICA



HIGH GRADE FUSED SILICA



Bundle Diameter (A) inches (mm)	Total Length (B) inches (mm)	Numerical Aperture	F/#	Acceptance Cone 2a (degrees)	Model No.
Glass Fiber Optic Bundles (400 to 1500 nm)					
0.063 (1.6)	48 (1219)	0.56	0.9	68	77522
0.125 (3.2)	24 (610)	0.56	0.9	68	77524
0.125 (3.2)	36 (914)	0.56	0.9	68	77525
0.125 (3.2)	48 (1219)	0.56	0.9	68	77526
0.125 (3.2)	72 (1829)	0.56	0.9	68	77527
Standard Grade Fused Silica Fiber Optic Bundles (280 to 2200 nm)					
0.125 (3.2)	24 (610)	0.27	1.9	31	77563
0.125 (3.2)	36 (914)	0.27	1.9	31	77564
High Grade Fused Silica Fiber Optic Bundles (260 to 2200 nm)					
0.063 (1.6)	24 (610)	0.22	2.3	25	77575
0.063 (1.6)	36 (914)	0.22	2.3	25	77576
0.125 (3.2)	24 (610)	0.22	2.3	25	77577