

Interface with: M-UMR5 and M-MVN50  
4 holes M3 clr / □ 20

Interface with: M-UTR46 and M-PO46  
4 holes M3 clr / Ø 38

Interface with: M-UMR8 and M-MVN80  
4 holes M4 clr / □ 34

2 holes M4 clr, 40 spacing

Interface: M-UTR80 with M-BR4  
4 holes M4 thd / □ 63

Interface with: M-UTR80 and M-PO80  
4 holes M4 clr / Ø 68 (or □ 48)

4 holes M4 thd / □ 108

2 holes M5 clr, 40 spacing

Interface with: M-UMR12 and M-MVN120  
4 holes M5 clr / □ 92

Interface with: M-UTR120  
4 holes M5 clr / Ø 105

Interface with: Rod and Optics table  
1 hole M6 clr

