



Model	Thread A	Dimensions [in. (mm)]	
		B	C
Male-Male Thread Adaptors			
TA-8Q20-10	8-32	1/4-20	0.25 (6.4)
TA-8M6-10	8-32	M6	0.25 (6.4)
TA-8M4-10	8-32	M4	0.21 (5.3)
TA-M4M6-10	M4	M6	0.25 (6.4)
TA-M4Q20-10	M4	1/4-20	0.25 (6.4)
TA-M6Q20-10	M6	1/4-20	0.21 (5.3)