









SECTION A-A



MODEL: LPXY-1H OR M-LPXY-1H



MODEL	D in (mm)	MOUNTING HOLE THREAD
LPXY-1H	0.50 (12.7)	8-32
M-LPXY-1H	0.47 (12.0)	M4 x 0.7



TECH SPECS:

1. X / Y TRANSLATION: $\pm 1.6\text{mm}$ AT $254\mu\text{m/rev}$; $0.71\mu\text{m}$ SENSITIVITY
2. CLEAR APERTURE: Ø.90 [22.9mm]
3. OPTICS AXIS HEIGHT: 1.0 [25.4mm]
4. MAX OPTICS EDGE THICKNESS: 0.40 [10.2mm]